



MSL 8 Digital

High Volume Magnetic Stirrer

Be Assured. Choose MSL Series



Safety, Resistance and Reliability.

- Clear digital display to accurately set up the desired working conditions
- Gentlestart and low speed stirring to protect your sample
- Brushless motor ensures premium power and long life span for up to 8L volume
- Resistant stainless steel cover allowing easy visualization of particles and color changes
- Constant stirring from 40 to 1500 rpm with SpeedServo technology
- MSL 8 Digital remains cold during continuous stirring
- Integraltimer and programmable auto-reverse of the stirring direction
- Safety Lock avoids possible unintentional changes
- Controlpanel is shielded from liquid spills
- The intermittent mode option enhances programming capabilities to match any application requirement



INSTRUMENT - CODE

MSL 8 Digital

100-240V/50-60Hz

F203A0510



TECHNICAL DATA

| | |
|---------------------------------|-------------------|
| USER INTERFACE | Digital |
| MOTOR | Brushless Motor |
| PLATE MATERIAL | Stainless steel |
| STIRRING VOLUME (H2O) | Up to 8 L |
| STIRRING SPEED RANGE | 40-1500 rpm |
| SPEED ACCURACY | 10 rpm |
| SPEED CONTROL | Digital |
| TORQUE COMPENSATION | SpeedServo™ |
| TIMER | Yes |
| AUTO-REVERSE | Yes |
| LOCK FUNCTION | Yes |
| INTERMITTENT MODE | Yes |
| AUTOMATIC SHUT DOWN | Yes, after 10 sec |
| DIMENSION (WXHxD) | 168x60x183 mm |
| PERMISSIBLE AMBIENT TEMPERATURE | 5-40°C |
| POWER INPUT | 5 W |
| WEIGHT | 0.95 Kg |



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
155 Keyland Court, Bohemia
NY 11716 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +86 18616509163
velpchina@velp.com